

An easy finish with long-term benefits.

Calibra Universal cement lets you achieve easy cleanup with enduring results.

75% of dentists have had to use a bur to remove excess resin cement¹



Achieve immediate and reliable bond strength

Calibra Universal cement offers easy cleanup and long-lasting results with no bonding agent required. A longer tack cure window of up to 10 seconds and 45-second extended gel phase gives clinicians the time needed for thorough excess cement cleanup.

- Self-adhesive – no separate bonding agent required
- High immediate 6 minute bond strength for a strong initial bond
- Can be used for nearly any indications including oxide ceramics (Zirconia), glass ceramics (lithium silicates), composite and metal-based inlays, onlays crowns, bridges, and posts
- 5 esthetic shade options with Shade Stable™ Technology virtually eliminates color shift over time
- High radiopacity

Robust Bonding Performance

Get the reliability and predictability of proven immediate and short term bond strength for long-lasting results.



Tips for Success

- Tooth preparation should leave the surface evenly glistening with moisture
- Ensure the tooth preparation is a retentive preparation. The axial angle should have a taper between 8°-10° and the cervical-occlusal height should be greater than 4mm.
- LED lights like SmartLite® Focus® with a single peak output of 470nm are recommended
- High power, dual, or broad spectrum lights may cause premature hardening of excess cement – check curing light effect on mixed cement in the laboratory prior to clinical use
- Remove floss horizontally through interproximals during cleanup so as not to dislodge the restoration before cement has completely set
- Following all excess removal, exposed margins may be light cured 20-40 seconds to assist restoration stabilization



1. The Key Group, Inc.: 2015 Dental Omnibus: Quarter 1 Voice of Customer Reports

2. 10-second tack cure window equals five-second wave cure per surface. For excess cement cleanup, monowave output LED lights with a single peak output around 470nm are recommended. High power, dual or broad spectrum lights may cause premature

RelyX Unicem 2 and Maxcem Elite are not registered trademarks of DENTSPLY.

Ordering Information

Dual Cure AutoMix Syringe Refill Package

2 Syringes (4.5g ea); 20 Mixing Tips

607402 Light Shade

607403 Medium Shade

607405 Translucent Shade

607406 Opaque Shade

607407 Bleach Shade

607086 Dual Cure AutoMix Syringe Mixing Tip Refill (50)

